

兔抗 HOXA11/D11 多克隆抗体

中文名称：兔抗 HOXA11/D11 多克隆抗体

英文名称：Anti-HOXA11/D11 rabbit polyclonal antibody

别名：HOX1; HOX11/HOX4; HOX4F

抗原：HOXA11/D11

储存：冷冻（-20℃）避光

宿主：Rabbit

反应种属：Human Mouse

相关类别：一抗

标记物：Unconjugate

克隆类型：Unconjugate

技术规格

Background:	Sequence-specific transcription factor which is part of a developmental regulatory system that provides cells with specific positional identities on the anterior-posterior axis.
Applications:	WB, IF
Name of antibody:	HOXA11/D11
Immunogen:	Synthesized peptide derived from internal of human HOXA11/D11.
Full name:	homeobox A11/D11
Synonyms :	HOX1; HOX11/HOX4; HOX4F
SwissProt:	P31270/P31277
WB Predicted band size:	34/35 kDa

WB Positive control:	RAW264.7 cells lysate
WB Recommended dilution :	500-3000
IF positive control:	Huvec cells
IF Recommend dilution:	100-500



