

HP 抗原（重组蛋白）

中文名称：HP 抗原（重组蛋白）

英文名称：HP Antigen (Recombinant Protein)

储 存：冷冻（-20℃）

相关类别：抗原

概述

Fusion protein corresponding to a region derived from 162-406 amino acids of human HP

技术规格

Full name:	haptoglobin
Synonyms:	BP, HPA1S, HP2ALPHA2
Swissprot:	P00738
Gene Accession:	BC121125
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a preproprotein, which is processed to yield both alpha and beta chains, which subsequently combine as a tetramer to produce haptoglobin. Haptoglobin functions to bind free plasma hemoglobin, which allows degradative enzymes to gain access to the hemoglobin, while at the same time preventing loss of iron through the kidneys and protecting the kidneys from damage by hemoglobin. Mutations in this gene and/or its regulatory regions cause ahaptoglobinemia or hypohaptoglobinemia. This gene has also been linked to diabetic nephropathy, the incidence of coronary artery disease in type 1 diabetes, Crohn's disease

e, inflammatory disease behavior, primary sclerosing cholangitis, susceptibility to idiopathic Parkinson's disease, and a reduced incidence of Plasmodium falciparum malaria. A similar duplicated gene is located next to this gene on chromosome 16. Multiple transcript variants encoding different isoforms have been found for this gene.