

## 兔抗 GULP1 多克隆抗体

中文名称：兔抗 GULP1 多克隆抗体

英文名称：Anti-GULP1 rabbit polyclonal antibody

别名：GULP, engulfment adaptor PTB domain containing 1; CED6; GULP; CED-6

相关类别：一抗

抗原：GULP1

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene is an adapter protein necessary for the engulfment of apoptotic cells by phagocytes. Several transcript variants, some protein coding and some thought not to be protein coding, have been found for this gene.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	GULP1
<b>Immunogen:</b>	Synthetic peptide of human GULP1

<b>Full name:</b>	GULP, engulfment adaptor PTB domain containing 1
<b>Synonyms:</b>	CED6; GULP; CED-6
<b>SwissProt:</b>	Q9UBP9
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human liver cancer and human esophagus cancer
<b>IHC Recommend dilution:</b>	50-

