

兔抗 GAS1 多克隆抗体

中文名称: 兔抗 GAS1 多克隆抗体

英文名称: Anti-GAS1 rabbit polyclonal antibody

储 存: 冷冻 (-20℃)

抗 原: GAS1

宿 主: Rabbit

反应种属: Human, Mouse

相关类别: 一抗

标 记 物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	Growth arrest-specific 1 plays a role in growth suppression. GAS1 blocks entry to S phase and prevents cycling of normal and transformed cells. Gas1 is a putative tumor suppressor gene.
Applications:	ELISA, WB
Name of antibody:	GAS1
Immunogen:	Synthetic peptide of human GAS1
Full name:	growth arrest-specific 1
SwissProt:	P54826
ELISA Recommended dilution:	1000-2000
WB Predicted band size:	36 kDa
WB Positive control:	Mouse brain tissue
WB Recommended dilution:	200-1000



