

## 兔抗 HOXB4 多克隆抗体

中文名称：兔抗 HOXB4 多克隆抗体

英文名称：Anti-HOXB4 rabbit polyclonal antibody

别名：HOX2; HOX2F; HOX-2.6

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：HOXB4

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene is a member of the Antp homeobox family and encodes a nuclear protein with a homeobox DNA-binding domain. It is included in a cluster of homeobox B genes located on chromosome 17. The encoded protein functions as a sequence-specific transcription factor that is involved in development. Intracellular or ectopic expression of this protein expands hematopoietic stem and progenitor cells in vivo and in vitro, making it a potential candidate for therapeutic stem cell expansion.
<b>Applications:</b>	ELISA, WB
<b>Name of antibody:</b>	HOXB4
<b>Immunogen:</b>	Synthetic peptide of human HOXB4

<b>Full name:</b>	homeobox B4
<b>Synonyms :</b>	HOX2; HOX2F; HOX-2.6
<b>SwissProt:</b>	P17483
<b>ELISA Recommended dilution:</b>	1000-2000
<b>WB Predicted band size:</b>	28 kDa
<b>WB Positive control:</b>	RAW264.7 cells
<b>WB Recommended dilution:</b>	200-1000

