

兔抗 GYS1 多克隆抗体

中文名称：兔抗 GYS1 多克隆抗体

英文名称：Anti-GYS1 rabbit polyclonal antibody

别名：GSY; GYS

相关类别：一抗

储存：冷冻（-20℃）

抗原：GYS1

宿主：Rabbit

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The protein encoded by this gene catalyzes the addition of glucose monomers to the growing glycogen molecule through the formation of alpha-1,4-glycoside linkages. Mutations in this gene are associated with muscle glycogen storage disease. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.
Applications:	ELISA, WB, IHC
Name of antibody:	GYS1

Immunogen:	Fusion protein of human GYS1
Full name:	glycogen synthase 1
Synonyms :	GSY; GYS
SwissProt:	P13807
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human brain and human thyroid cancer
IHC Recommend dilution:	50-200
WB Predicted band size:	84 kDa
WB Positive control:	293T and Hela cells
WB Recommended dilution:	200-1000



