

兔抗 GYS2 多克隆抗体

中文名称：兔抗 GYS2 多克隆抗体

英文名称：Anti-GYS2 rabbit polyclonal antibody

相关类别：一抗

抗 原：GYS2

储 存：冷冻（-20℃）

宿 主：Rabbit

反应种属：Human, Mouse, Rat

标 记 物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The protein encoded by this gene, liver glycogen synthase, catalyzes the rate-limiting step in the synthesis of glycogen - the transfer of a glucose molecule from UDP-glucose to a terminal branch of the glycogen molecule. Mutations in this gene cause glycogen storage disease type 0 (GSD-0) - a rare type of early childhood fasting hypoglycemia with decreased liver glycogen content.
Applications:	ELISA, IHC
Name of antibody:	GYS2
Immunogen:	Fusion protein of human GYS2

Full name:	glycogen synthase 2
SwissProt:	P54840
ELISA Recommended dilution:	2000-10000
IHC positive control:	Human ovarian cancer and human brain
IHC Recommend dilution:	100-300

