

兔抗 GAS2L3 多克隆抗体

中文名称：兔抗 GAS2L3 多克隆抗体

英文名称：Anti-GAS2L3 rabbit polyclonal antibody

别名：G2L3

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：GAS2L3

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Gas2L3 (growth arrest-specific 2 like 3) is a 694 amino acid protein that contains one calponin-homology (CH) domain and may function in a similar manner to Gas2. The gene encoding Gas2L3 maps to human chromosome 12, which encodes over 1,100 genes and comprises approximately 4.5% of the human genome. Cytoskeletal linker protein. May promote and stabilize the formation of the actin and microtubule network.
Applications:	ELISA, IHC
Name of antibody:	GAS2L3
Immunogen:	Synthetic peptide of human GAS2L3

Full name:	growth arrest specific 2 like 3
Synonyms:	G2L3
SwissProt:	Q86XJ1
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human cervical cancer
IHC Recommend dilution:	10-50

