

## 兔抗 HOXB6 多克隆抗体

中文名称：兔抗 HOXB6 多克隆抗体

英文名称：Anti-HOXB6 rabbit polyclonal antibody

别名：HOX2; HU-2; HOX2B; Hox-2.2

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：HOXB6

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：Unconjugate

### 技术规格

<b>Background:</b>	This gene is a member of the Antp homeobox family and encode s a protein with a homeobox DNA-binding domain. It is included in a cluster of homeobox B genes located on chromosome 17. The encoded protein functions as a sequence-specific transcription factor that is involved in development, including that of lung and skin, and has been localized to both the nucleus and cytoplasm. Altered expression of this gene or a change in the subcellular localization of its protein is associated with some cases of acute myeloid leukemia and colorectal cancer.
<b>Applications:</b>	ELISA, IHC

<b>Name of antibody:</b>	HOXB6
<b>Immunogen:</b>	Synthetic peptide of human HOXB6
<b>Full name:</b>	homeobox B6
<b>Synonyms :</b>	HOX2; HU-2; HOX2B; Hox-2.2
<b>SwissProt:</b>	P17509
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human brain and human gastric cancer



