

IRAK3 抗原（重组蛋白）

中文名称：IRAK3 抗原（重组蛋白）

英文名称：IRAK3 Antigen (Recombinant Protein)

别名：ASRT5; IRAKM

储存：冷冻（-20℃）

相关类别：抗原

概述

Fusion protein corresponding to a region derived from 165-452 amino acids of human IRAK3

技术规格

Full name:	interleukin-1 receptor-associated kinase 3
Synonyms:	ASRT5; IRAKM
Swissprot:	Q9Y616
Gene Accession:	BC057800
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a member of the interleukin-1 receptor-associated kinase protein family. Members of this family are essential components of the Toll/IL-R immune signal transduction pathways. This protein is primarily expressed in monocytes and macrophages and functions as a negative regulator of Toll-like receptor signaling. Mutations in this gene are associated with a susceptibility to asthma. Alternate splicing results in multiple transcript variants.