

KRT4 抗原（重组蛋白）

中文名称：KRT4 抗原（重组蛋白）

英文名称：KRT4 Antigen (Recombinant Protein)

储 存： 冷冻（-20℃）

相关类别： 抗原

概述

Fusion protein corresponding to a region derived from 151-461 amino acids of human KRT4

技术规格

| | |
|---------------------------|--|
| Full name: | keratin 4, type II |
| Synonyms: | K4; CK4; CK-4; CYK4; WSN1 |
| Swissprot: | P19013 |
| Gene Accession: | BC042174 |
| Purity: | >85%, as determined by Coomassie blue stained SDS-PAGE |
| Expression system: | Escherichia coli |
| Tags: | His tag C-Terminus, GST tag N-Terminus |
| Background: | The protein encoded by this gene is a member of the keratin gene family. The type II cytokeratins consist of basic or neutral proteins which are arranged in pairs of heterotypic keratin chains coexpressed during differentiation of simple and stratified epithelial tissues. This type II cytokeratin is specifically expressed in differentiated layers of the mucosal and esophageal epithelia with family member KRT13. |