

兔抗 GBP4 多克隆抗体

中文名称：兔抗 GBP4 多克隆抗体

英文名称：Anti-GBP4 rabbit polyclonal antibody

别名：guanylate binding protein 4; Mpa2

相关类别：一抗

抗原：GBP4

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Guanylate-binding proteins, such as GBP4, are induced by interferon and hydrolyze GTP to both GDP and GMP (Vestal, 2005 [PubMed 16108726]). Plays a role in erythroid differentiation.
Applications:	ELISA, IHC
Name of antibody:	GBP4
Immunogen:	Synthetic peptide of human GBP4
Full name:	guanylate binding protein 4
Synonyms:	Mpa2
SwissProt:	Q96PP9
ELISA Recommended diluti	5000-10000

on:	
IHC positive control:	Human cervical cancer and Human colorectal cancer
IHC Recommend dilution:	30-150



