

## 兔抗 HASPIN 多克隆抗体

中文名称：兔抗 HASPIN 多克隆抗体

英文名称：Anti-HASPIN rabbit polyclonal antibody

别名：histone H3 associated protein kinase; GSG2

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：HASPIN

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	Serine/threonine-protein kinase that phosphorylates histone H3 at 'Thr-3' (H3T3ph) during mitosis. This positions and activates AURKB and other components of the chromosomal passenger complex (CPC) at centromeres to ensure proper chromatid cohesion, metaphase alignment and normal progression through the cell cycle.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	HASPIN
<b>Immunogen:</b>	Synthetic peptide of human HASPIN

<b>Full name:</b>	histone H3 associated protein kinase
<b>Synonyms:</b>	GSG2
<b>SwissProt:</b>	Q8TF76
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human tonsil and Human liver cancer
<b>IHC Recommend dilution:</b>	40-200



