

## SNCG 抗原（重组蛋白）

中文名称：SNCG 抗原（重组蛋白）

英文名称： SNCG Antigen (Recombinant Protein)

储 存： 冷冻（-20℃）

相关类别： 抗原

概述

Fusion protein corresponding to a region derived from 71-451 amino acids of human SNCG

技术规格

<b>Full name:</b>	synuclein, gamma (breast cancer-specific protein 1)
<b>Synonyms:</b>	SR; BCSG1
<b>Swissprot:</b>	O76070
<b>Gene Accession:</b>	BC014098
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	This gene encodes a member of the synuclein family of proteins which are believed to be involved in the pathogenesis of neurodegenerative diseases. Mutations in this gene have also been associated with breast tumor development.