

兔抗 HAUS2 多克隆抗体

中文名称：兔抗 HAUS2 多克隆抗体

英文名称：Anti-HAUS2 rabbit polyclonal antibody

别名：HAUS augmin like complex subunit 2; CEP27; C15orf25; HsT17025

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：HAUS2

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a subunit of the augmin complex. The augmin complex plays a role in microtubule attachment to the kinetochore and central spindle formation.
Applications:	ELISA, WB
Name of antibody:	HAUS2
Immunogen:	Synthetic peptide of human HAUS2
Full name:	HAUS augmin like complex subunit 2
Synonyms:	CEP27; C15orf25; HsT17025

SwissProt:	Q9NVX0
ELISA Recommended dilution:	5000-10000
WB Predicted band size:	27 kDa
WB Positive control:	Mouse lung tissue lysate
WB Recommended dilution:	200-1000

