

TBX22 抗原（重组蛋白）

中文名称：TBX22 抗原（重组蛋白）

英文名称：TBX22 Antigen (Recombinant Protein)

储 存： 冷冻（-20℃）

相关类别： 抗原

概述

Fusion protein corresponding to a region derived from 321-520 amino acids of human TBX22

技术规格

Full name:	T-box 22
Synonyms:	CPX, CLPA, TBXX, dJ795G23.1
Swissprot:	Q9Y458
Gene Accession:	BC014194
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene is a member of a phylogenetically conserved family of genes that share a common DNA-binding domain, the T-box. T-box genes encode transcription factors involved in the regulation of developmental processes. Mutations in this gene have been associated with the inherited X-linked disorder, Cleft palate with ankyloglossia, and it is believed to play a major role in human palatogenesis. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.