

TNFRSF12A 抗原（重组蛋白）

中文名称：TNFRSF12A 抗原（重组蛋白）

英文名称：TNFRSF12A Antigen (Recombinant Protein)

储 存：冷冻（-20℃）

相关类别：抗原

概述

Fusion protein corresponding to a region derived from 28-129 amino acids of human TNFRSF12A

技术规格

Full name:	Tumor necrosis factor receptor superfamily, member 12A
Synonyms:	FN14; CD266; TWEAKR
Swissprot:	Q9NP84
Gene Accession:	BC002718
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	Tumor necrosis factor receptor superfamily member 12A also known as the Tweak-receptor (TweakR) is a protein that in humans is encoded by the TNFRSF12A gene. Receptor for TNFSF12/TWEAK. Weak inducer of apoptosis in some cell types. Promotes angiogenesis and the proliferation of endothelial cells. May modulate cellular adhesion to matrix proteins. Highly expressed in heart, placenta and kidney. Intermediate expression in lung, skeletal muscle and pancreas.