

兔抗 HSD17B14 多克隆抗体

中文名称：兔抗 HSD17B14 多克隆抗体

英文名称：Anti-HSD17B14 rabbit polyclonal antibody

别名：DHRS10; SDR47C1; retSDR3

相关类别：一抗

抗原：HSD17B14

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	17-beta-hydroxysteroid dehydrogenases, such as HSD17B14, are primarily involved in metabolism of steroids at the C17 position and also of other substrates, such as fatty acids, prostaglandins, and xenobiotics.
Applications:	ELISA, WB, IHC
Name of antibody:	HSD17B14
Immunogen:	Fusion protein of human HSD17B14
Full name:	hydroxysteroid (17-beta) dehydrogenase 14
Synonyms :	DHRS10; SDR47C1; retSDR3
SwissProt:	Q9BPX1
ELISA Recommended dilution:	1000-2000

IHC positive control:	Human colon cancer and human liver cancer
IHC Recommend dilution:	50-200
WB Predicted band size:	28 kDa
WB Positive control:	Mouse kidney tissue
WB Recommended dilution:	200-1000



