

兔抗 MYBPC1 多克隆抗体

中文名称: 兔抗 MYBPC1 多克隆抗体

英文名称: Anti-MYBPC1 rabbit polyclonal antibody

别名: LCCS4; MYBPCC; MYBPCS

储存: 冷冻 (-20℃) 避光

抗原: MYBPC1

宿主: Rabbit

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: Unconjugate

技术规格

Background:

This gene encodes a member of the myosin-binding protein C family. Myosin-binding protein C family members are myosin-associated proteins found in the cross-bridge-bearing zone (C region) of A bands in striated muscle. The encoded protein is the slow skeletal muscle isoform of myosin-binding protein C and plays an important role in muscle contraction by recruiting muscle-type creatine kinase to myosin filaments. Mutations in this gene are associated with distal arthrogryposis type I. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.

Applications:	WB
Name of antibody:	MYBPC1
Immunogen:	Synthesized peptide derived from internal of human MYBP C1.
Full name:	myosin binding protein C, slow type
Synonyms :	LCCS4; MYBPCC; MYBPCCS
SwissProt:	Q00872
WB Predicted band size:	128 kDa
WB Positive control:	293 cells lysate
WB Recommended dilution:	500-3000

