

ENO2 抗原（多肽片段）

中文名称： ENO2 抗原 （多肽片段）

英文名称： ENO2 Antigen

相关类别： 抗原

储 存： 冷冻（-20℃）

概述

Synthetic peptide corresponding to a region derived from 420-434 amino acids of human enolase 2 (gamma, neuronal)

技术规格

Full name	Enolase 2 (gamma, neuronal)
Synonyms	NSE
Swissprot	P09104
Gene Accession	NP_001966
Background	This gene encodes one of the three enolase isoenzymes found in mammals. This isoenzyme, a homodimer, is found in mature neurons and cells of neuronal origin. A switch from alpha enolase to gamma enolase occurs in neural tissue during development in rats and primates. Has neurotrophic and neuroprotective properties on a broad spectrum of central nervous system (CNS) neurons. Binds, in a calcium-dependent manner, to cultured neocortical neurons and promotes cell survival.