

ERGIC3 抗原（重组蛋白）

中文名称： ERGIC3 抗原（重组蛋白）

英文名称： ERGIC3 Antigen (Recombinant Protein)

别名： ERGIC and golgi 3; Erv46; CGI-54; C2orf47; PRO0989; C20orf47; NY-BR-84; SDBCAG84; dJ477O4.2

储存： 冷冻（-20℃）

相关类别： 抗原

概述

Fusion protein corresponding to a region derived from 184-383 amino acids of human ERGIC3

技术规格

Full name:	ERGIC and golgi 3
Synonyms:	Erv46; CGI-54; C2orf47; PRO0989; C20orf47; NY-BR-84; SDBCAG84; dJ477O4.2
Swissprot:	Q9Y282
Gene Accession:	BC014014
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	Possible role in transport between endoplasmic reticulum and Golgi