

兔抗 HSPG2 多克隆抗体

中文名称： 兔抗 HSPG2 多克隆抗体

英文名称： Anti-HSPG2 rabbit polyclonal antibody

别名： PLC; SJA; SJS; HSPG; SJS1; PRCAN

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： HSPG2

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes the perlecan protein, which consists of a core protein to which three long chains of glycosaminoglycans (heparan sulfate or chondroitin sulfate) are attached. The perlecan protein is a large multidomain proteoglycan that binds to and cross-links many extracellular matrix components and cell-surface molecules. It has been shown that this protein interacts with laminin, prolargin, collagen type IV, FGFBP1, FBLN2, FGF7 and transthyretin, etc., and it plays essential roles in multiple biological activities.
Applications:	ELISA, IHC

Name of antibody:	HSPG2
Immunogen:	Synthetic peptide of human HSPG2
Full name:	heparan sulfate proteoglycan 2
Synonyms :	PLC; SJA; SJS; HSPG; SJS1; PRCAN
SwissProt:	P98160
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human thyroid cancer and Human brain
IHC Recommend dilution:	50-200

