

兔抗 NAA15 多克隆抗体

中文名称：兔抗 NAA15 多克隆抗体

英文名称：Anti-NAA15 rabbit polyclonal antibody

别名：Ga19; NATH; TBDN; NARG1; NAT1P; TBDN100

抗原：NAA15

储存：冷冻（-20℃）避光

宿主：Rabbit

反应种属：Human Mouse

相关类别：一抗

标记物：Unconjugate

克隆类型：Unconjugate

技术规格

Background:	This gene encodes a protein of unknown function. However , similarity to proteins in yeast and other species suggests that this protein may be an N-acetyltransferase.
Applications:	WB
Name of antibody:	NAA15
Immunogen:	Synthesized peptide derived from internal of human NARG1 .
Full name:	N(alpha)-acetyltransferase 15, NatA auxiliary subunit
Synonyms :	Ga19; NATH; TBDN; NARG1; NAT1P; TBDN100
SwissProt:	Q9BXJ9
WB Predicted band size:	101 kDa

WB Positive control:	Huvec cells lysate
WB Recommended dilution:	500-3000

