

## 兔抗 PYGM 多克隆抗体

中文名称：兔抗 PYGM 多克隆抗体

英文名称：Anti-PYGM rabbit polyclonal antibody

相关类别：一抗

储 存：冷冻（-20℃）

宿 主：Rabbit

抗 原：PYGM

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a muscle enzyme involved in glycogenolysis. Highly similar enzymes encoded by different genes are found in liver and brain. Mutations in this gene are associated with McArdle disease (myophosphorylase deficiency), a glycogen storage disease of muscle. Alternative splicing results in multiple transcript variants.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	PYGM
<b>Immunogen:</b>	Fusion protein of human PYGM
<b>Full name:</b>	phosphorylase, glycogen, muscle
<b>SwissProt:</b>	P11217

<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human thyroid cancer and Human liver cancer
<b>IHC Recommend dilution:</b>	30-150
<b>WB Predicted band size:</b>	97 kDa
<b>WB Positive control:</b>	Mouse heart tissue
<b>WB Recommended dilution:</b>	500-2000



