

## 兔抗 RAB6A 多克隆抗体

中文名称：兔抗 RAB6A 多克隆抗体

英文名称：Anti-RAB6A rabbit polyclonal antibody

别名：RAB6

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：RAB6A

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

#### Background:

This gene encodes a member of the RAB family, which belongs to the small GTPase superfamily. GTPases of the RAB family bind to various effectors to regulate the targeting and fusion of transport carriers to acceptor compartments. This protein is located at the Golgi apparatus, which regulates trafficking in both a retrograde (from early endosomes and Golgi to the endoplasmic reticulum) and an anterograde (from the Golgi to the plasma membrane) directions. Myosin II is

	<p>s an effector of this protein in these processes. This protein is also involved in assembly of human cytomegalovirus (HCMV) by interacting with the cellular protein Bicaudal D1, which interacts with the HCMV virion tegument protein, pp150. Multiple alternatively spliced transcript variants encoding different isoforms have been identified.</p>
<b>Applications:</b>	ELISA, WB
<b>Name of antibody:</b>	RAB6A
<b>Immunogen:</b>	Synthetic peptide of human RAB6A
<b>Full name:</b>	RAB6A, member RAS oncogene family
<b>Synonyms :</b>	RAB6
<b>SwissProt:</b>	P20340
<b>ELISA Recommended dilution:</b>	1000-2000
<b>WB Predicted band size:</b>	24 kDa
<b>WB Positive control:</b>	293T, A549 and A172 cells
<b>WB Recommended dilution:</b>	200-1000

