

兔抗 TRO 多克隆抗体

中文名称：兔抗 TRO 多克隆抗体

英文名称：Anti-TRO rabbit polyclonal antibody

别名：trophinin; MAGED3; MAGE-d3

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：TRO

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

| | |
|--------------------------|--|
| Background: | This gene encodes a membrane protein that mediates cell adhesion between trophoblastic cells and the epithelial cells of the endometrium. The encoded protein participates in cell signalling during embryo implantation, and may also be involved in cancer formation. This gene is located near several other closely related genes on chromosome X. Alternative splicing results in multiple transcript variants. |
| Applications: | ELISA, IHC |
| Name of antibody: | TRO |
| Immunogen: | Fusion protein of human TRO |
| Full name: | trophinin |

| | |
|------------------------------------|-----------------|
| Synonyms: | MAGED3; MAGE-d3 |
| SwissProt: | Q12816 |
| ELISA Recommended dilution: | 5000-10000 |
| IHC positive control: | Human tonsil |
| IHC Recommend dilution: | 50-100 |

