

兔抗 RAC3 多克隆抗体

中文名称：兔抗 RAC3 多克隆抗体

英文名称：Anti-RAC3 rabbit polyclonal antibody

别名：ras-related C3 botulinum toxin substrate 3 (rho family, small GTP binding protein Rac3)

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：RAC3

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a GTPase which belongs to the RAS superfamily of small GTP-binding proteins. Members of this superfamily appear to regulate a diverse array of cellular events, including the control of cell growth, cytoskeletal reorganization, and the activation of protein kinases. Alternative splicing results in multiple transcript variants.
Applications:	ELISA, WB, IHC
Name of antibody:	RAC3
Immunogen:	Full length fusion protein

Full name:	ras-related C3 botulinum toxin substrate 3 (rho family, small GTP binding protein Rac3)
SwissProt:	P60763
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil
IHC Recommend dilution:	50-300
WB Predicted band size:	21 kDa
WB Positive control:	Human empty ileal tissue and Human ileum tissue lysates
WB Recommended dilution:	500-2000



