

兔抗 RALGAPA1 多克隆抗体

中文名称：兔抗 RALGAPA1 多克隆抗体

英文名称：Anti-RALGAPA1 rabbit polyclonal antibody

别名：Ral GTPase activating protein catalytic alpha subunit 1; p240; GRIPE; GARNL1; TULIP1; RalGAPalpha1

抗原：RALGAPA1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a major subunit of the RAL-GTPase activating protein. A similar protein in mouse binds E12, a transcriptional regulator of immunoglobulin genes. The mouse protein also functions in skeletal muscle by binding to the regulatory 14-3-3 proteins upon stimulation with insulin or muscle contraction. A pseudogene of this gene has been identified on chromosome 9.
Applications:	ELISA, IHC
Name of antibody:	RALGAPA1

Immunogen:	Synthetic peptide of human RALGAPA1
Full name:	Ral GTPase activating protein catalytic alpha sub unit 1
Synonyms:	p240; GRIPE; GARNL1; TULIP1; RalGAPalpha1
SwissProt:	Q6GYQ0
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human cervical cancer and Human thyroid cancer
IHC Recommend dilution:	40-200



