

兔抗 RAPGEF1 多克隆抗体

中文名称: 兔抗 RAPGEF1 多克隆抗体

英文名称: Anti-RAPGEF1 rabbit polyclonal antibody

别名: Rap guanine nucleotide exchange factor 1; C3G; GRF2

抗原: RAPGEF1

储存: 冷冻 (-20℃)

宿主: Rabbit

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	This gene encodes a human guanine nucleotide exchange factor. It transduces signals from CRK by binding the SH3 domain of CRK, and activating several members of the Ras family of GTPases. This signaling cascade that may be involved in apoptosis, integrin-mediated signal transduction, and cell transformation. Several alternatively spliced transcript variants of this gene have been described, but the full-length nature of some variants has not been determined.
Applications:	ELISA, IHC
Name of antibody:	RAPGEF1

Immunogen:	Synthetic peptide of human RAPGEF1
Full name:	Rap guanine nucleotide exchange factor 1
Synonyms:	C3G; GRF2
SwissProt:	Q13905
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil
IHC Recommend dilution:	25-100



