

小鼠抗 mBanana 单克隆抗体

中文名称： 小鼠抗 mBanana 单克隆抗体

英文名称： Anti-mBanana mouse monoclonal antibody

相关类别： 标签抗体

储 存： 冷冻（-20℃） 避光

宿 主： Mouse

抗 原： mBanana

反应种属： mBanana

标 记 物： Unconjugate

克隆类型： mouse monoclonal

技术规格

Background:	mBanana is a novel monomeric red fluorescent protein mutant. mBanana is a Fruit Fluorescent Protein which was developed in Dr. Roger Tsien' s lab by directed mutagenesis of mRFP1, a monomeric mutant of DsRed . Fruit Fluorescent Proteins such as mBanana have demonstrated stable expression and perform successfully in numerous fusion applications.
Applications:	WB
Name of antibody:	mBanana
Immunogen:	Purified mBanana protein expressed in HEK293 cell
WB Predicted band size:	27 kDa
WB Positive control:	HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY mBanana



	ne) cDNA for 48 hrs and lysed
WB Recommended dilution:	1000-5000