

兔抗 RAPSN 多克隆抗体

中文名称：兔抗 RAPSN 多克隆抗体

英文名称：Anti-RAPSN rabbit polyclonal antibody

别名：receptor associated protein of the synapse; FADS; CMS11; CMS4C; RAPSYN; RNF205

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：RAPSN

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:

This gene encodes a member of a family of proteins that are receptor associated proteins of the synapse. The encoded protein contains a conserved cAMP-dependent protein kinase phosphorylation site, and plays a critical role in clustering and anchoring nicotinic acetylcholine receptors at synaptic sites by linking the receptors to the underlying postsynaptic cytoskeleton, possibly by direct association with actin or spectrin. Mutations in this gene may play a role in postsynaptic

	congenital myasthenic syndromes. Alternatively s pliced transcript variants encoding multiple isofo rms have been observed for this gene.
Applications:	ELISA, IHC
Name of antibody:	RAPSN
Immunogen:	Fusion protein of human RAPSN
Full name:	receptor associated protein of the synapse
Synonyms:	FADS; CMS11; CMS4C; RAPSYN; RNF205
SwissProt:	Q13702
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	25-50

