

兔抗 ACTB 多克隆抗体

中文名称：兔抗 ACTB 多克隆抗体

英文名称： Anti-ACTB rabbit polyclonal antibody

别名： BRWS1; PS1TP5BP1

相关类别： 内参抗体

储存： 冷冻（-20℃）

抗原： ACTB

宿主： Rabbit

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes one of six different actin proteins. Actins are highly conserved proteins that are involved in cell motility, structure, and integrity. This actin is a major constituent of the contractile apparatus and one of the two nonmuscle cytoskeletal actins.
Applications:	ELISA, WB, IHC
Name of antibody:	ACTB
Immunogen:	Synthetic peptide of human ACTB
Full name:	actin, beta
Synonyms:	BRWS1; PS1TP5BP1

SwissProt:	P60709
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human lung cancer and Human tonsil
IHC Recommend dilution:	50-200
WB Predicted band size:	42 kDa
WB Positive control:	Mouse brain tissue, Rat brain tissue, Rat liver tissue, Mouse liver tissue, Human cerebrum tissue lysates
WB Recommended dilution:	1000-5000



