

兔抗 NDUFB10 多克隆抗体

中文名称：兔抗 NDUFB10 多克隆抗体

英文名称： Anti-NDUFB10 rabbit polyclonal antibody

别 名： PDSW

相关类别： 一抗

储 存： 冷冻（-20℃）

宿 主： Rabbit

抗 原： NDUFB10

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	NADH dehydrogenase [ubiquinone] 1 beta subcomplex subunit 10 is an enzyme that in humans is encoded by theNDUFB10 gene. NADH dehydrogenase (ubiquinone) 1 beta subcomplex subunit 10 is an accessory subunit of the NADH dehydrogenase (ubiquinone) complex, located in the mitochondrial inner membrane. It is also known as Complex I and is the largest of the five complexes of the electron transport chain.
Applications:	ELISA, WB, IHC
Name of antibody:	NDUFB10
Immunogen:	Full length fusion protein

Full name:	NADH dehydrogenase (ubiquinone) 1 beta subcomplex, 10, 22kDa
Synonyms :	PDSW
SwissProt:	O96000
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human breast cancer and human esophagus cancer
IHC Recommend dilution:	25-100
WB Predicted band size:	21 kDa
WB Positive control:	Jurkat and A549 cell, human fetal muscle tissue
WB Recommended dilution:	1000-5000



