

## 兔抗 ACTIN(植物)多克隆抗体

中文名称：兔抗 ACTIN(植物)多克隆抗体

英文名称：Anti-ACTIN(Plants) rabbit polyclonal antibody

相关类别：内参抗体

储 存：冷冻 (-20℃)

宿 主：Rabbit

抗 原：ACTIN

反应种属：Plants

标 记 物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	Member of actin subclass composed of ACT12 and ACT4. RNA is expressed at very low levels in vegetative organs, low levels in flowers and very high levels in pollen. Expression of an ACT12/GUS fusion was found in vascular tissues, tapetum, developing and mature pollen, the root cap and in a ring of pericycle tissues during lateral root initiation and early development.
<b>Applications:</b>	ELISA, WB
<b>Name of antibody:</b>	ACTIN
<b>Immunogen:</b>	Synthetic peptide of Arabidopsis thaliana actin
<b>Full name:</b>	Actin
<b>SwissProt:</b>	P53497

<b>ELISA Recommended dilution:</b>	5000-10000
<b>WB Predicted band size:</b>	42 kDa
<b>WB Positive control:</b>	Nicotiana tabacum leaves lysate and Corn leaves lysates
<b>WB Recommended dilution:</b>	500-2000

