

兔抗 RBBP5 多克隆抗体

中文名称：兔抗 RBBP5 多克隆抗体

英文名称：Anti-RBBP5 rabbit polyclonal antibody

别名：RBQ3; SWD1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：RBBP5

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a ubiquitously expressed nuclear protein which belongs to a highly conserved subfamily of WD-repeat proteins. The encoded protein binds directly to retinoblastoma protein, which regulates cell proliferation. It interacts preferentially with the underphosphorylated retinoblastoma protein via the E1A-binding pocket B. Three alternatively spliced transcript variants that encode different protein isoforms have been described for this gene.
Applications:	ELISA, WB
Name of antibody:	RBBP5

Immunogen:	Fusion protein of human RBBP5
Full name:	retinoblastoma binding protein 5
Synonyms :	RBQ3; SWD1
SwissProt:	Q15291
ELISA Recommended dilution:	1000-2000
WB Predicted band size:	59 kDa
WB Positive control:	Hela, hepg2 and HT-29 cells
WB Recommended dilution:	200-1000

