

兔抗 IDS 多克隆抗体

中文名称： 兔抗 IDS 多克隆抗体

英文名称： Anti-IDS rabbit polyclonal antibody

别名： iduronate 2-sulfatase; MPS2; SIDS

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： IDS

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a member of the sulfatase family of proteins. The encoded preproprotein is proteolytically processed to generate two polypeptide chains. This enzyme is involved in the lysosomal degradation of heparan sulfate and dermatan sulfate. Mutations in this gene are associated with the X-linked lysosomal storage disease mucopolysaccharidosis type II, also known as Hunter syndrome. Alternative splicing results in multiple transcript variants, at least one of which encodes a preproprotein that is proteolytically processed.
Applications:	ELISA, IHC
Name of antibody:	IDS

Immunogen:	Fusion protein of human IDS
Full name:	iduronate 2-sulfatase
Synonyms:	MPS2; SIDS
SwissProt:	P22304
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	50-200

