

## 兔抗 TUSC1 多克隆抗体

中文名称：兔抗 TUSC1 多克隆抗体

英文名称：Anti-TUSC1 rabbit polyclonal antibody

别名：TSG9; TSG-9

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：TUSC1

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene is located within the region of chromosome 9p that harbors tumor suppressor genes critical in carcinogenesis. It is an intronless gene which is downregulated in non-small-cell lung cancer and small-cell lung cancer cell lines, suggesting that it may play a role in lung tumorigenesis.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	TUSC1
<b>Immunogen:</b>	Synthetic peptide of human TUSC1
<b>Full name:</b>	tumor suppressor candidate 1

<b>Synonyms :</b>	TSG9; TSG-9
<b>SwissProt:</b>	Q2TAM9
<b>ELISA Recommended dilution:</b>	1000-2000
<b>IHC positive control:</b>	Human liver cancer
<b>IHC Recommend dilution:</b>	30-150
<b>WB Predicted band size:</b>	23 kDa
<b>WB Positive control:</b>	PC3 cell and human colon cancer tissue
<b>WB Recommended dilution:</b>	200-1000



