

AP 标记兔抗小鼠 IgG

中文名称： AP 标记兔抗小鼠 IgG

英文名称： AP-conjugated Rabbit anti-mouse IgG

相关类别： AP 标记二抗

储 存： 冷冻（-20℃）

标 记 物： AP

技术规格

Background	Antibody Specificity: Based on immunoelectrophoresis and/or ELISA, the antibody reacts with whole molecule mouse IgG. It also reacts with the light chains of other mouse immunoglobulins. No antibody was detected against non-immunoglobulin serum proteins. The antibody may cross-react with immunoglobulins from other species.
Applications	1/2000-20000 for ELISA 1/2000-20000 for Western blotting