

Biotin 标记兔抗大鼠 IgG

中文名称: Biotin 标记兔抗大鼠 IgG

英文名称: Biotin-conjugated Rabbit anti-rat IgG

别 名: 生物素标记兔抗大鼠 IgG

相关类别: Biotin 标记二抗

储 存: 冷冻(-20℃)

技术规格

0

Background

Immunoglobulin G (IgG) is one of the most abundant proteins in serum with normal levels between 8-17 mg/mL in adult blo od. IgG is important for our defense against microorganisms a nd the molecules are produced by B lymphocytes as a part of our adaptive immune response. The IgG molecule has two sep arate functions, to bind to the pathogen that elicited the response and to recruit other cells and molecules to destroy the ant igen. The variability of the IgG pool is generated by somatic recombination and the number of specificities in an individual at a given time point is estimated to be 1011 variants.



Applications	1/1000-5000 for ELISA 1/1000-5000 for Western blotting 1/100-
	1000 for Immunohistochemistry/cytochemistry