

生物素标记山羊抗人 IgG

中文名称： 生物素标记山羊抗人 IgG

英文名称： Biotin-conjugated Goat anti-human IgG

相关类别： Biotin 标记二抗

储 存： 冷冻（-20℃）

标记物： Biotin

技术规格

Background	Antibody Specificity: Based on immunoelectrophoresis and/or ELISA , the antibody reacts with whole molecule mouse IgG. It also reacts with the light chains of other mouse immunoglobulins. No antibody was detected against non-immunoglobulin serum proteins. The antibody may cross-react with immunoglobulins from other species .
Applications	1/5000-100000 for ELISA and Western blotting using enzyme-conjugated streptavidin 1/100-1000 for enzyme immunohisto/cytochemis

try 1/50-200 for Flow Cytometer and fluorescence immunohisto/cyt
ochemistry