

## 生物素标记山羊抗兔 IgG

中文名称：生物素标记山羊抗兔 IgG

英文名称：Biotin-conjugated Goat anti-Rabbit IgG

相关类别：Biotin 标记二抗

储 存： 冷冻（-20℃）

标记 物： Biotin

### 技术规格

<b>Background</b>	Antibody Specificity: Based on immunoelectrophoresis and/or ELISA, the antibody reacts with whole molecule rabbit IgG. It also reacts with the light chains of other rabbit immunoglobulins. No antibody was detected against non-immunoglobulin serum proteins. The antibody may cross-react with immunoglobulins from other species.
<b>Applications</b>	1/2000-20000 for ELISA 1/2000-20000 for Western blotting