

## 生物素标记山羊抗小鼠 IgG

中文名称: 生物素标记山羊抗小鼠 IgG

英文名称: Biotin-conjugated Goat anti-mouse IgG

相关类别: Biotin 标记二抗

储 存: 冷冻(-20℃)

标记物: Biotin

## 技术规格

Background	Antibody Specificity: Based on immunoelectrophoresis and/or
	ELISA, the antibody reacts with whole molecule mouse IgG. It
	also reacts with the light chains of other mouse immunoglobu
	lins. No antibody was detected against non-immunoglobulin s
	erum proteins. The antibody may cross-react with immunoglob
	ulins from other species.
Applications	1/5000-100000 for ELISA and Western blotting using enzyme-c
	onjugated streptavidin 1/100-1000 for enzyme immunohisto/cy
	tochemistry 1/50-200 for Flow Cytometer and fluorescence im
	munohisto/cytochemistry

