

兔抗 TYK2 多克隆抗体

中文名称：兔抗 TYK2 多克隆抗体

英文名称：Anti-TYK2 rabbit polyclonal antibody

别名：tyrosine kinase 2; JTK1; IMD35

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：TYK2

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:

This gene encodes a member of the tyrosine kinase and, more specifically, the Janus kinases (JAKs) protein families. This protein associates with the cytoplasmic domain of type I and type II cytokine receptors and promulgate cytokine signals by phosphorylating receptor subunits. It is also component of both the type I and type III interferon signaling pathways. As such, it may play a role in anti-viral immunity. A mutation in this gene has been associated with hyperimmunoglobulin E syndrome (HIES) - a primary immunodeficiency characterized by elevated serum i

	mmunoglobulin E.
Applications:	ELISA, IHC
Name of antibody:	TYK2
Immunogen:	Fusion protein of human TYK2
Full name:	tyrosine kinase 2
Synonyms:	JTK1; IMD35
SwissProt:	P29597
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer and Human thyroid cancer
IHC Recommend dilution:	150-300



