

HRP 标记的兔抗鸡 lgG

- 中文名称: HRP标记的兔抗鸡 IgG
- 英文名称: HRP-conjugated Rabbit anti-chicken IgG
- 相关类别: HRP标记二抗
- 储 存: 冷冻 (-20℃)
- 标 记 物: HRP

技术规格

| Background | Immunoglobulin G (IgG), is one of the most abundant proteins in serum with normal levels between 9.5-12.5 mg/ml in adult blood. I gG is important for our defence against microorganisms and the m olecules are produced by B lymphocytes as a part of our adaptive immune response. The IgG molecule has two separate functions; to bind to the pathogen that elicited the response and to recruit othe r cells and molecules to destroy the antigen. The variability of the IgG pool is generated by Somatic recombination and the number o f specificities in an individual at a given time point is estimated to be 1011 variants. |
|--------------|---|
| Applications | 1/1000-5000 for ELISA 1/1000-5000 for Western blotting 1/200-500 for Immunohistochemistry/cytochemistry |