

## HRP 标记的小鼠抗兔 IgG，轻链特异性

中文名称：HRP 标记的小鼠抗兔 IgG，轻链特异性

英文名称：HRP-conjugated IgG Fraction monoclonal Mouse Anti-Rabbit IgG, Light chain specific

相关类别：HRP 标记二抗

储 存：冷冻（-20℃）

标 记 物：HRP

### 技术规格

<b>Background</b>	Based on immunoelectrophoresis and/or ELISA, the antibody reacts with whole molecule rabbit IgG. It also reacts with the light chains of other rabbit immunoglobulins. No antibody was detected against non-immunoglobulin serum proteins. The antibody has been tested by ELISA and/or solid-phase adsorbed to ensure minimal cross-reaction with bovine, chicken, goat, guinea pig, Armenian hamster, horse, human, mouse, rat, and sheep serum proteins, but it may cross-react with immunoglobulins from other species.
<b>Applications</b>	1/ 2000-50000 for ELISA 1/ 200-2000 for Immunohisto/cytochemistry 1/ 5000-100000 for Western blotting with ECL substrates