

HRP 标记山羊抗大鼠 IgG

中文名称： HRP 标记山羊抗大鼠 IgG

英文名称： HRP-conjugated Goat anti-rat IgG

相关类别： HRP 标记二抗

储 存： 冷冻 (-20℃)

标 记 物： HRP

技术规格

| | |
|---------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Background | Antibody Specificity: Based on immunoelectrophoresis and/or ELISA, the antibody reacts with whole molecule rat IgG. It also reacts with the light chains of other rat immunoglobulins. No antibody was detected against non-immunoglobulin serum proteins. The antibody may cross-react with immunoglobulins from other species. |
| Applications | 1/2000-50000 for ELISA 1/200-2000 for Immunohisto/cytochemistry 1/5000-100000 for Western blotting with ECL substrates |