

兔抗 IGF1R(Ab-1346) 多克隆抗体

中文名称：兔抗 IGF1R(Ab-1346) 多克隆抗体

英文名称：Anti-IGF1R(Ab-1346) rabbit polyclonal antibody

别名：IGFR; CD221; IGFIR; JTK13

抗原：IGF1R(Ab-1346)

储存：冷冻（-20℃）避光

宿主：Rabbit

反应种属：Human

相关类别：一抗

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This receptor binds insulin-like growth factor with a high affinity. It has tyrosine kinase activity. The insulin-like growth factor I receptor plays a critical role in transformation events. Cleavage of the precursor generates alpha and beta subunits. It is highly overexpressed in most malignant tissues where it functions as an anti-apoptotic agent by enhancing cell survival.
Applications:	WB, IHC
Name of antibody:	IGF1R(Ab-1346)
Immunogen:	Synthesized non-phosphopeptide derived from human IGF1R

	around the phosphorylation site of tyrosine 1346 (Q-P-Y(p)-A-H).
Full name:	insulin-like growth factor 1 receptor
Synonyms :	IGFR; CD221; IGFR; JTK13
SwissProt:	P08069
IHC positive control:	Human brain tissue
IHC Recommend dilution:	50-100
WB Predicted band size:	155 kDa
WB Positive control:	A549 cells lysate
WB Recommended dilution:	500-3000



