

兔抗 RELA(Phospho-Ser281) 多克隆抗体

中文名称：兔抗 RELA(Phospho-Ser281) 多克隆抗体

英文名称：Anti-RELA(Phospho-Ser281) rabbit polyclonal antibody

别名：p65; NFKB3

抗原：RELA(Phospho-Ser281)

储存：冷冻（-20℃）避光

宿主：Rabbit

相关类别：一抗

反应种属：Human Mouse Rat

标记物：Unconjugate

克隆类型：Unconjugate

技术规格

WB Recommended dilution:	500-1000
WB Positive control:	Mouse brain tissue lysates
WB Predicted band size:	60 kDa
SwissProt:	Q04206
Synonyms :	p65; NFKB3
Full name:	v-rel avian reticuloendotheliosis viral oncogene homolog A
Immunogen:	Peptide sequence around phosphorylation site of serine281(E-L-S(p)-E-P) derived from Human NF- kappa B .
Name of antibody:	RELA(Phospho-Ser281)

Applications:	WB
Back ground:	<p>NF-kappa-B is a ubiquitous transcription factor involved in several biological processes. It is held in the cytoplasm in an inactive state by specific inhibitors. Upon degradation of the inhibitor, NF-kappa-B moves to the nucleus and activates transcription of specific genes. NF-kappa-B is composed of NFKB1 or NFKB2 bound to either REL, RELA, or RELB. The most abundant form of NF-kappa-B is NFKB1 complexed with the product of this gene, RELA. Four transcript variants encoding different isoforms have been found for this gene.</p>

