

AMCA 标记的驴抗兔 IgG

中文名称： AMCA 标记的驴抗兔 IgG

英文名称： AMCA-conjugated Donkey Anti-Rabbit IgG

相关类别： 其他二抗

储 存： 冷冻 (-20℃)

标 记 物： AMCA

技术规格

Background	Based on immunoelectrophoresis and/or ELISA, the antibody reacts with whole molecule rabbit IgG. It also reacts with the light chains of other rabbit immunoglobulins. No antibody was detected against non- immunoglobulin serum proteins. The antibody has been tested by ELISA and/or solid-phase adsorbed to ensure minimal cross-reaction with bovine, chicken, goat, guinea pig, Syrian hamster, horse, human, mouse, rat, and sheep serum proteins, but it may cross-react with immunoglobulins from other species.
Applications	1/50-200 for most application