

兔抗 UBASH3A 多克隆抗体

中文名称: 兔抗 UBASH3A 多克隆抗体

英文名称: Anti-UBASH3A rabbit polyclonal antibody

别名: TULA; CLIP4; STS-2; TULA-1

抗原: UBASH3A

储存: 冷冻 (-20℃)

宿主: Rabbit

反应种属: Human, Mouse

标记物: Unconjugate

相关类别: 一抗

克隆类型: rabbit polyclonal

技术规格

Background:	This gene encodes one of two family members belonging to the T-cell ubiquitin ligand (TULA) family. Both family members can negatively regulate T-cell signaling. This family member can facilitate growth factor withdrawal-induced apoptosis in T cells, which may occur via its interaction with AIF, an apoptosis-inducing factor. Alternative splicing of this gene results in multiple transcript variants.
Applications:	ELISA, IHC
Name of antibody:	UBASH3A
Immunogen:	Fusion protein of human UBASH3A

Full name:	ubiquitin associated and SH3 domain containing A
Synonyms:	TULA; CLIP4; STS-2; TULA-1
SwissProt:	P57075
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil and human esophagus cancer
IHC Recommend dilution:	40-250



